

China Huawei World Huawei

Huawei adheres to the essence of Chinese traditional medical culture, combines with the introduction of new technology and innovation and creates health care products with contemporary significance, based on China and have the whole world in view.



Wuhan HuaWei Technology Co., Ltd was founded in 2002, it is a high-tech enterprise devoted to medical health products' development, production and sales, company adhering to "quality first, credibility supreme" as our business philosophy, Take "to serve the community, to benefit the human health" as our responsibility; After 10 years of development, company's main products are at the advanced level in the industry, and have established good relations of cooperation and Stable cooperation relationship with the domestic and foreign customers.



EXCELLENT · CONTINUED · SHARING

Every Detail Of Life,
We Guard



Wuhan Huawei Technology Co.,Ltd.

Office Address: 8/F, Building A, Spaceflight Star Town, Two Seven Road, Jiang'an District, 430012, Wuhan, Hubei, China

Factory Address: B-F11-3, Huazhong International Industrial Park, Yangluo Port, 430415, Wuhan, Hubei, China

Tel: +86-27-88868920 **Fax:** +86-27-88876095 **Website:** www.whhuawei.com.cn

Email: jimmywang@whhuawei.com.cn sales3@whhuawei.com.cn


sales4@whhuawei.com.cn

 **Wuhan Huawei Technology Co.,Ltd.**

Company Profile



Workshop views



Company located beside the scenic Yangtze River, has a modern production plant in Huazhong International Industrial Park, Wuhan Yangluo Port, where is a high-tech enterprise agglomeration with modern and advanced medical materials production line, PU composite material production line, external plaster production line, lotion production line, cream and other pharmaceutical production lines, also packing equipment, company is the research base of Chinese Medical Association of Chinese Medicine for "curing internal disease from external" and the high-tech products' production base of the Wuhan University Medical School, in addition to the production of company's own proprietary brand, it becomes more and more high-end products' OEM base .


Wuhan Huawei Technology Co., Ltd. was founded in 2002, it is a high-tech enterprise devoted to medical health products' development, production and sales, company adheres to "quality first, credibility supreme" as our business philosophy, takes "to serve the community, to benefit the human health" as our responsibility, after 10 years of development, company's main products are at the advanced level in the industry, and has established good relations and stable cooperation with domestic and foreign customers .

Company has an own independent and modernized R&D center, it formed a joint open development system with universities in Beijing, Shanghai, Wuhan, Zhengzhou and other related R&D institutions, it has maintained close cooperation with Switzerland FORBO, U.S. Rohm Haas' R&D experts as well. After years of research and development, company has successes to develop the Huawei breathable membrane (waterborne polyurethane imitation skin membrane, we called it WPU), it is domestic initiative product. Its water resistance, anti-bacteria, breathable, no sensitization properties which are completely beyond the traditional film material, It is the world's leading technical level and has obtained the national invention patent .

Glance of certificates



Win-win cooperation



Wuhan Huawei Technology Co., Ltd. is committed to external using products, we have created a new upgraded production process, quality and efficacy, our products embody the wisdom and efforts of Swiss scientist, professional team of medical plaster industry, Rohm and Haas senior experts, and the experts from Wuhan University who has studied in US, we use the world first human body effect method to test the transdermal effect of the patch plaster. Our infrared patch series, medical dressing series, breathable film series, special ointment series, gynecological gel series are exported to Europe and the United States, Russia, Southeast Asia, have gained customers' trust and praises from both domestic and overseas .

For customers at home and abroad, we offer the following services.

1. Wholesales, agents and distribution business: anyone who has the sales strength and channel, if interest in our products, we can grant the distribution agents to reach a win-win cooperation .
2. Any customer, as long as the volume of orders meet our MOQ, we can provide OEM, ODM service, spare you from the tedious work of designing, production and other aspects of work, you can focus on market development with all your energy and time .

We have been guided by laws of the market. Consistent keep the same pace with our distributors, and are responsible for all the patients, do our best for national health, and had gained strong support and recognition from all parties. We have many successful operation of the market available for your reference. We will maintain a pragmatic attitude and philosophy of win-win cooperation .

Contents

目 录

Advanced Wound Care Series	01
Silicone Foam Wound Dressing	03
Hydrocolloid Wound Dressing	06
Hydrogel Wound Dressing	09
Alginate Wound Dressing	10
Wound Care Series	11
Surgical Incise Dressing	12
I.V. Cannula Fixing Dressing	18
Transparent Wound Dressing	19
Nonwoven Wound Dressing	21
Adhesive Bandage and Fixing Series	24
Catheter Fixing Dressing	25
Medical Tapes	28
Silicone Adhesive Tape	31
Wound Skin Closure Strip	32
Personal Healthcare Series	34
Silicone Scar Gel Sheet	35
Acne Patch	37
Pain Relief Patches	39
Cold Therapy Products	42
Functional Herbal Patches	44
Personal Protective Equipment	48
Isolation Gowns	49
Face Shield	50
Disposable Face Mask	51
Medical Raw Material Jumbo Roll Series	52
PU Material	54
PE Material	55
Nonwoven and Woven Material	56

WUHAN HUAWEI TECHNOLOGY CO.,LTD.


Advanced Wound Care Series

Silicone Foam Wound Dressing
Hydrocolloid Wound Dressing
Hydrogel Wound Dressing
Alginate Wound Dressing

CE FDA ISO 9001 ISO 13485

SILICONE FOAM WOUND DRESSING

Silicone foam dressing is an absorbent, atraumatic dressing made from absorbent core and silicone gel. The wound contact surface of Silicone foam dressing is coated with a layer of soft silicone gel that does not stick to the surface of a wound or cause trauma to delicate new tissue upon removal. It has follow advantage:
1, Wet Healing . 2, Breathable and skin-friendly . 3, Soft to check the wound . 4, Reduce pain . 5, Can be pasted repeatedly .



Indications

The product is suitable for many types of exuding wounds including leg and pressure ulcers superficial and partial thickness burns, donor sites, postoperative wounds and skin abrasions. The dressing absorbs exudate and maintains a moist wound-healing environment while minimising the risk of maceration .

Characteristics

Moist Healing:

The three-layer structure provides a dynamic exudate and wet wound healing environment that promotes rapid wound healing.

Breathable and skin-friendly:

The porous structure facilitates the conduction of the osmotic fluid and has a large absorption capacity for the liquid and rapidly absorbs the effusion.

Soft to check the wound:

The material is soft, comfortable, non-stickable to hair, and does not damage newborn epithelium and granulation. The wound contact layer can be reattached after peeling off the dressing to observe the wound without losing its adhesiveness.

Reduce pain:

Suitable for sensitive populations, children, cancer patients and chronic unhealed wounds. Mild tack can reduce wound damage and pain caused by changing dressings.

Imported raw materials:

This product is made of imported high-absorbency foam, silicone gel and separating film

Specification



Silicone wound dressing (thinner boarder type) .



Foam Dressing (Laminated with PU, with silicone adhesive at back) .

Item No.	Size(cm)	Packing
11.01.09	18*18	10pcs/box,24boxes/CTN
11.01.15	14*14	10pcs/box,24boxes/CTN

Item No.	Size(cm)	Packing
11.02.01	10*10	10pcs/box,36boxes/CTN
11.02.02	15*15	10pcs/box,24boxes/CTN

WUHAN HUAWEI TECHNOLOGY CO.,LTD.


SILICONE FOAM WOUND DRESSING





Silicone wound dressing (thinner boarder type) .

Item No.	Size(cm)	Packing
11.01.03	10*15	10pcs/box,36boxes/CTN
11.01.04	10*20	10pcs/box,36boxes/CTN
11.01.05	10*25	10pcs/box,30boxes/CTN
11.01.06	10*30	10pcs/box,30boxes/CTN



Silicone wound dressing (thinner boarder type) .

Item No.	Size(cm)	Packing
11.01.01	7.5*7.5	10pcs/box,54boxes/CTN
11.01.02	10*10	10pcs/box,36boxes/CTN
11.01.07	12.5*12.5	10pcs/box,36boxes/CTN
11.01.08	15*15	10pcs/box,24boxes/CTN
11.01.12	20*20	10pcs/box,40boxes/CTN



Silicone wound dressing (thinner boarder type) .

Item No.	Size(cm)	Packing
11.01.11	18.5*24	10pcs/box,20boxes/CTN



Silicone wound dressing (thinner boarder type) .

Item No.	Size(cm)	Packing
11.01.10	16*20	10pcs/box,20boxes/CTN



Foam Dressing (Laminated with PU, Non adhesive).

Item No.	Size(cm)	Packing
11.03.01	10*10	10pcs/box,36boxes/CTN
11.03.02	15*15	10pcs/box,24boxes/CTN



Tracheostomy Foam Dressing

Item No.	Size(cm)	Packing
11.04.01	10*10	10pcs/box,36boxes/CTN
11.04.04	7*10	10pcs/box,36boxes/CTN



WUHAN HUAWEI TECHNOLOGY CO.,LTD.

HYDROCOLLOID WOUND DRESSING



HYDROCOLLOID WOUND DRESSING

Hydrocolloid dressing is a kind of elastic dressing made of medical hot-melt adhesive with CMC(carboxymethylcellulose) hydrophilic granules.



Product description

1.The semi-permeable matte polyurethane or polyurethane foam to enhance the absorbency.Film thickness:0.025-0.035mm.

2.When the thickness of hydrocolloid daessing is 0.4mm thick:

Moisture permeability: 0.34g/10cm²/24h/37°C

1.06g/10cm²/72h/37°C

Hygroscopic capacity: 2.74g/10cm²/24h/37°C

3.52g/10cm²/72h/37°C(only for reference)

3.Hygroscopic capacity: 2.83g/g/24h/37°C

3.23g/g/72h/37°C(only for reference)

4.PH 5.7 Acid & Alkali Resistance PH 3-12

5.Hygroscopic capacity and hygroscopic capability can be abjusted according needs.

Hygroscopic ability: 2.40g/10cm²/24h/37

2.46g/10cm²/72h/37°C

Product performance

1.Excellent exudates absorption capability.

2.Keep the wound in a moist environment, promote wound Healing,effectively release pains and reduce wound care frequency.

3.Water-proof protect the wound from microbial invasion.

4.Excellent elasticith and self-adhesive,convenient and comfortable to wear.

5.Please change the dressing when it turns white.

6.Easy to peel off, no adherence to wound, so effectively reduce the potential damage and pains.

7.Available in different sizes and shapes to meet the clinical needs.

Application

For slight or moderate exudating wounds

1. Leg ulcer. 2. I - II degree bedsores 3. Superficial and burn wound. 4. Post-operation wound.

5. Skin-donor sites. 6. All kinds of superficial wound and cosmetic wound.

7. For granulation and epithelialisation process of chronic wound.

Cautions

1. Make sure the wound area is clean and dry before using dressing.

2. Hydrocolloid dressing should be 2 cm bigger than the wound area to make sure that daessing could be cover the wound.

3. If the depth of wound is more than 5 mm, it is better to fill the wound with proper material before using hydrocolloid dressing.

4. It is not suitable for wounds with heavy exudates.

5. When the dressing becomes white and swelling, it indicates that dressing should be changed.


6. At the very beginning of using the dressing, wound area maybe enlarged, this is caused by the debridement function of the dressing and it is normal phenomenon.

7. The gel will be formed by the mixture of hydrocolloid molecule and exudates. As it looks like purulence secretion, it would be misunderstood as the wound infection. Just clean it with saline.

8. There maybe some smell from the dressing sometimes, and this smell disappear after cleaning the wound with saline water.

9. Dressing should be changed at once when leakage occurred around the dressing.

Specifications




Ultra thin

Item No.	Size(cm)	Packing	CTN(cm)
04.18.11	5*5	20pcs/box, 20boxes/CTN	43*17*19
04.18.13	10*10	10pcs/box, 20boxes/CTN	22*22*24
04.18.14	15*15	10pcs/box, 20boxes/CTN	21*21*17
04.18.15	5*10	10pcs/box, 20boxes/CTN	21*21*17
04.18.16	5*15	10pcs/box, 20boxes/CTN	21*21*17



With thin border

Item No.	Size(cm)	Packing	CTN(cm)
04.18.04	5*5	20pcs/box, 20boxes/CTN	27.5*21.5*21
04.18.06	10*10	10pcs/box, 20boxes/CTN	43*17*19
04.18.07	15*15	10pcs/box, 20boxes/CTN	43*17*19



With thin border
(For foot healing)

Item No.	Size(cm)	Packing	CTN(cm)
04.18.08	8*12	10pcs/box, 20boxes/CTN	43*17*19




With thin border
(For sacrum)

Item No.	Size(cm)	Packing	CTN(cm)
04.18.09	15*18	10pcs/box, 20boxes/CTN	22*22*24



With thin border

Item No.	Size(cm)	Packing	CTN(cm)
04.18.12	4.2*6.2	6pcs/box, 120boxes/CTN	43*17*19
04.18.17	4.4*6.9	6pcs/box, 120boxes/CTN	43*17*19



With thin border
(For sacrum)

Item No.	Size(cm)	Packing	CTN(cm)
04.18.10	14*18	10pcs/box, 20boxes/CTN	22*22*24



With thin border

Item No.	Size(cm)	Packing	CTN(cm)
04.18.21	4.2*6.2	6pcs/box, 120boxes/CTN	43*17*19



With thin border

Item No.	Size(cm)	Packing	CTN(cm)
04.18.18	2*6	6pcs/box, 120boxes/CTN	43*17*19
04.18.19	2.5*6.5	6pcs/box, 120boxes/CTN	43*17*19
04.18.20	6*1.7	6pcs/box, 120boxes/CTN	43*17*19

HYDROGEL WOUND DRESSING

Hydrogel wound dressing (also called aquogel wound dressing) is a new type of wound dressing consisting of a high water-content absorbent polymer gel .


Performance

1. Debridement gel contains more than 90% water and it has the dual function of moisturizing the wound and absorbing the seepage which can keep the wound moist and effectively avoid the dry necrosis of the wound tissue .
2. Debridement gel plays the exogenous debridement effect to wrap, separate and clean the wound necrotic tissue, bacteria and so on .
3. Improving the regeneration ability of wound granulation tissue, promoting epithelial cell division and migration , accelerating the healing of the wound, and avoiding to leave a scar .
4. It's soft, comfortable and elastic which can reduce pain and avoid adhesive and secondary damage to the wound .
5. The cool sensation provided by the hydrogel to the wound offers relief from pain for at least six hours .

Application

1. Dry or dehydrated wounds
2. Partial or full-thickness lesions
3. Abrasions or severe scrapes
4. Minor burns, I - II degree burned wound and skin donor site.
5. Wounds with granulated tissue development
6. Radiation skin damage

Specifications



Item No.	Size(cm)	Packing	CTN(cm)	Remarks
10.01.01	5*5	10pcs/box, 54boxes/CTN	45.5*39.5*33	Hydrogel wound dressing(self adhesive)
10.01.02	10*10	10pcs/box, 36boxes/CTN	45.5*33*35	
10.01.03	15*15	10pcs/box, 24boxes/CTN	43*41*35	



Item No.	Size(cm)	Packing	CTN(cm)	Remarks
10.02.01	7.5*10	10pcs/box, 20boxes/CTN	43*29*32.5	Hydrogel wound dressing(with adhesive boarder)
10.02.02	10*15	10pcs/box, 20boxes/CTN	42.5*29*33	
10.02.03	10*20	10pcs/box, 20boxes/CTN	53*29*33	
10.02.04	10*25	10pcs/box, 20boxes/CTN	57*31.5*32.5	
10.02.05	12.5*12.5	10pcs/box, 20boxes/CTN	51.5*41*35	

WUHAN HUAWEI TECHNOLOGY CO.,LTD.

HYDROGEL WOUND DRESSING



ALGINATE WOUND DRESSING

Alginate dressing is the medical dressing of fiber and ion from the natural seaweed .

Product description

1. Good absorbency and gelation capacity .
2. Viscosity coefficient>800cp, Asb<0.5%, fiber denier:1.5D, strength:1.5g/D, Elongation>10%.
3. Water absorption≥ 10, bacteriostasis ratio>80%. Pieces and strips are available .

Characteristics

- 1.It is the mixture of natural seaweed and calcium ion, and it has good biocompatibility .
- 2.Gel could be formed when the dressing meets exudates and blood, it could protect the would area and promote the would healing .
- 3.The dressing has excellent absorbency, it is soft and comfortable .
- 4.Release of Calcium ion could activate prothrombin, and could stop bleeding effectively .
- 5.No adherence to the wound, it will protect the nerve terminals, reduce pain and easy to be peeled off .
- 6.Maceration of the skin around will not occur .
- 7.It is biologically degradable and environment-friendly .
- 8.It is soft and can be used to fill the wound, and enhance the growth of the cavity .

Application

- 1.All kinds of moderate and heavy exudated wounds .
- 2.All kinds of acute and chronic wounds hemorrhagic wounds .
- 3.All kinds of difficult healing wounds such as leg ulcer, bed sore, diabetic ulcer, post-operative wounds of tumour and abscess .
- 4.Skin-donor sites .
- 5.All kinds of lacuna, antrnm wounds, post-operative wounds of cavity sinus, exodontia .


Cautions

1. It is suitable for dry wound .
2. Make sure the wound area is clean and dry before using dressing .
3. Alginate dressing should be 2cm bigger than the wound area to make sure that dressing could cover skin completely .
4. The dressing should be stayed on wound less than one week .
5. When the exudates decrease, it suggests changing to another kind of dressing, such as hydrocolloid dressing .
6. Make sure the size, depth of the cavity wound before using, fill the wound from bottom. Don't left dead space, or may affect the wound healing .

Specifications



Item No.	Size(cm)	Packing	Remarks
12.01.01	5*5	20pcs/box 20boxes/CTN	120g Alginate dressing
12.01.02	10*10	10pcs/box 20boxes/CTN	120g Alginate dressing
12.01.03	2*30		120g Alginate fill dressing



WUHAN HUAWEI TECHNOLOGY CO.,LTD.

Wound Care Series

Surgical Incise Dressing

I. V. Cannula Dressing


Transparent Wound Dressing

Nonwoven Wound Dressing

CE FDA ISO 9001 ISO 13485

SURGICAL INCISE DRESSING

Medical surgical incise dressing film from Wuhan Huawei Technology co., Ltd. is divided into PU surgical film, PE surgical film, PU surgical iodine film, PE surgical iodine film. It adopt the high permeable polyurethane film (PU) or polyethylene (PE) as raw material, coated with a special medical pressure sensitive adhesive (iodine available), and then attached to release protective paper, protect the surgical area avoid cross-infection .



Structure and Characteristics

Hypoallergenic: Reduce the rate of sensitization to prevent infection .

Good Viscous: good adhesion, high-quality medical pressure-sensitive adhesion makes an good using effect without warped edge, no residue, no pain when you tear off, can be tightly retained for a long time on the skin .

Breathable and skin-friendly: highly breathable polyurethane film transparent and breathable to prevent moisture accumulation under surgery, provide a sterile working environment around the surgical incision .

Waterproof & resistance bacteria: iodine concentration of surgical iodine drape 0.45% -0.8%, effectively kill a variety of bacteria, sterilization rate reach 99%.

Humanized design: Easy-pull bars on both sides are available according to customers' requirements, or leave empty side (no glue) to tear off easier. It can also be used as the cut film with non-woven surgical drapes for compounding use, multiple sizes available .



Specifications

PE INCISE DRESSING



Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.14.01-E	15*15	50pcs/box 40boxes/CTN	60*22*33	Gamma-ray
04.14.02-E	10*20	30pcs/box 40boxes/CTN	60*22*33	Gamma-ray
04.14.03-E	20*30		58*28*30	Gamma-ray
04.14.04-E	30*30	20pcs/box 20boxes/CTN	73*28*30	Gamma-ray
04.14.05-E	30*40		73*28*30	Gamma-ray
04.14.06-E	30*45		73*28*30	Gamma-ray
04.14.07-E	45*45	20pcs/box 10boxes/CTN	58*28*30	Gamma-ray

PE IODINE INCISE DRESSING



Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.14.01-D	15*20	30pcs/box 40boxes/CTN	60*22*33	Gamma-ray
04.14.02-D	20*30	30pcs/box 20boxes/CTN	58*28*30	Gamma-ray
04.14.03-D	30*40	20pcs/box 20boxes/CTN	73*28*30	Gamma-ray
04.14.04-D	30*45		73*28*30	Gamma-ray

WUHAN HUAWEI TECHNOLOGY CO.,LTD.
SURGICAL INCISIVE DRESSING




PU IODINE INCISE DRESSING



Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.14.01-I	15*15	50pcs/box 40boxes/CTN	60*22*33	Gamma-ray
04.14.02-I	10*20	30pcs/box 20boxes/CTN	58*28*30	Gamma-ray
04.14.03-I	20*30	25pcs/box 20boxes/CTN	58*28*30	Gamma-ray
04.14.04-I	25*30		58*38*30	Gamma-ray
04.14.05-I	30*40		73*28*30	Gamma-ray
04.14.06-I	30*45	20pcs/box 10boxes/CTN	73*28*30	Gamma-ray
04.14.07-I	45*45		58*28*30	Gamma-ray
04.14.08-I	40*50		58*28*30	Gamma-ray
04.14.09-I	40*60	20pcs/box 20boxes/CTN	58*28*30	Gamma-ray
04.14.10-I	45*60		58*28*30	Gamma-ray
04.14.11-I	15*35	20pcs/box 20boxes/CTN	60*22*33	Gamma-ray
04.14.12-I	7*11	50pcs/box 20boxes/CTN	43*29*32.5	Gamma-ray

PU INCISE DRESSING



Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.14.01	15*15	50pcs/box 40boxes/CTN	60*22*33	Gamma-ray
04.14.02	10*20	30pcs/box 20boxes/CTN	58*28*30	Gamma-ray
04.14.03	20*30	25pcs/box 20boxes/CTN	58*28*30	Gamma-ray
04.14.04	25*30		58*38*30	Gamma-ray
04.14.05	30*40		73*28*30	Gamma-ray
04.14.06	30*45	20pcs/box 10boxes/CTN	73*28*30	Gamma-ray
04.14.07	45*45		58*28*30	Gamma-ray
04.14.08	40*50		58*28*30	Gamma-ray
04.14.09	40*60	20pcs/box 20boxes/CTN	58*28*30	Gamma-ray
04.14.10	45*60		58*28*30	Gamma-ray
04.14.11	15*35	20pcs/box 20boxes/CTN	60*22*33	Gamma-ray
04.14.12	7*11	50pcs/box 20boxes/CTN	43*29*32.5	Gamma-ray

SURGICAL INCISE DRESSING 3M TYPE WITH IODINE



Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.14.01-M	15*20	50pcs/box 40boxes/CTN	60*22*33	Gamma-ray
04.14.02-M	34*35	30pcs/box 40boxes/CTN	60*22*33	Gamma-ray
04.14.03-M	56*60	20pcs/box 20boxes/CTN	73*28*30	Gamma-ray
04.14.04-M	56*45		73*28*30	Gamma-ray
04.14.05-M	56*85	20pcs/box 10boxes/CTN	73*48*30	Gamma-ray
04.14.06-M	60*45	20pcs/box 20boxes/CTN	73*28*30	Gamma-ray
04.14.07-M	60*85	20pcs/box 10boxes/CTN	73*48*30	Gamma-ray

INCISE DRESSING WITH SINGLE WASTAGE COLLECTION BAG




Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.15.01S	15*20	20pcs/box 20boxes/CTN	53*32*33	Gamma-ray
04.15.02S	14*35		53*32*33	Gamma-ray
04.15.03S	30*40	10pcs/box 20boxes/CTN	60*36.5*30.5	Gamma-ray
04.15.04S	30*45		60*36.5*30.5	Gamma-ray
04.15.05S	45*45	10pcs/box 10boxes/CTN	53*32*33	Gamma-ray

DOUBLE WASTAGE COLLECTION BAG




Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.15.03D	30*40	10pcs/box 20boxes/CTN	60*36.5*30.5	Gamma-ray
04.15.04D	30*45		60*36.5*30.5	Gamma-ray
04.15.05D	45*45	10pcs/box 10boxes/CTN	53*32*33	Gamma-ray

THREE SIDE CONNECTION WASTAGE COLLECTION BAG



Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.15.05N	45*45	10pcs/box 10boxes/CTN	53.5*38*44.5	Gamma-ray
04.15.05I (with iodine)	45*45		53.5*38*44.5	Gamma-ray

FOUR SIDE CONNECTION WASTAGE COLLECTION BAG



Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.15.04N	30*40	10pcs/box 10boxes/CTN	53.5*38*44.5	Gamma-ray
04.15.04I (with iodine)	30*40		53.5*38*44.5	Gamma-ray

OEM SURGICAL DRAPES WITH INCISE DRESSING



Remarks: 1.All above shape is selectable our standard size. Customized size and shape are welcomed .
2.We can do based on customer size, shape, package, imprint, OEM order welcomed .


I.V. CANNULA FIXING DRESSING

Medical transparent wound dressing adopts high permeable aqueous polyurethane glue, processes through a special coating process, breathe together with skin .

Advantages

Green environmental protection, super breathable, waterproof anti-bacteria. Hypoallergenic, comfortable and convenient, single-permeability .

Polyurethane film
+ Hypoallergenic adhesive



Specifications



For fixing I.V. basic type
(with U port) .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.03.03	6*7	100pcs/box, 36boxes/CTN	53*32*23.5



for fixing carotid I.V catheter
(basic type) .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.03.04	6*7	100pcs/box, 36boxes/CTN	53*32*23.5
04.03.05	9*10	100pcs/box, 24boxes/CTN	53*32*23.5
04.03.06	10*15	100pcs/box, 20boxes/CTN	53*32*23.5
04.03.07	10*20	100pcs/box, 20boxes/CTN	53*32*23.5
04.03.08	10*25	100pcs/box, 10boxes/CTN	53*32*23.5



For fixing PICC and eye
dressings .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.03.11	6*7	50pcs/box, 40boxes/CTN	40.5*21*25.6



For fixing Y or butterfly
catheter or catheter
with drug kettle .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.03.10	6*7	50pcs/box 20boxes/CTN	40.5*21*25.5



For fixing PICC(basic) .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.03.09	10*12	100pcs/box, 24boxes/CTN	43*31*34



For fixing arterial and
T-extended catheter

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.03.12	7*8.5	50pcs/box, 20boxes/CTN	22*18.7*12.5



For fixing multi-lumen
catheter, pulmonary artery
catheter, hemodialysis
catheter, jujuar jujualar
vein and carotid catheter .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.03.13	8.5*10.5	50pcs/box, 24boxes/CTN	22*15*26.5



For donor sites and
anesthetic catheter fixing .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.03.16	10*15	100pcs/box, 24boxes/CTN	52*31*30
04.03.17	10*20	100pcs/box, 24boxes/CTN	58*31*30
04.03.18	10*25	100pcs/box, 24boxes/CTN	35*31*30
04.03.19	15*20	100pcs/box, 24boxes/CTN	41*30.5*29

WUHAN HUAWEI TECHNOLOGY CO.,LTD.

I.V. CANNULA FIXING DRESSING



TRANSPARENT WOUND DRESSING

Medical transparent dressing adopts contemporary international high-tech technique of new ultra-thin comfortable dressing, which is a special semipermeable membrane, with good breathability, hypoallergenic, high elasticity and other characteristics.



Structure and Characteristics

Waterproof & resistance bacteria: outer waterproof, breathable film to prevent moisture and bacteria invasion, reduce wound infection, good waterproof ensure that patients can bathe and other basic life activities.

Skin-friendly breathable: Adopt ultra-thin PU film to make the dressing with ultra-high moisture permeability, avoid moisture aggregation between the skin and the PU membrane, prolong use time, make the skin breathe freely, bestly care the injured.

High elastic PU film: Available to stick to any part of the body, reducing surrounding skin tension, comfortable and natural.

Super absorbent: absorbent pad with good absorb performance and greatly reduce the incidence of impregnation, provide a good moist healing environment to the wound.

Hypoallergenic: does not irritate the wound, good fixation, does not damage the skin when exposing.

Safe and easy to tear: the fluid guide net surface of absorbent pad can make the dressing not to stick with the wound, easy to peel when changing dressing, do not cause secondary damage to the wound.

Humanized design: complete specifications, more product styles, we can do according to customers' requirements, provide more choices for clinical use.

Instructions

1. Clean the wound, select the appropriate specifications of the dressing, remove the sterile packaging and peel off release paper.
2. Align absorbent pad to wound, gently affix it to the affected area then thoroughly strip release paper.
3. Remove the transparent protective film above, leave only dressing film.
4. Paste the transparent dressing around smoothly.



Application Postoperative wound, acute and chronic injury wound, small cut wound and scratches, etc.

Clinical departments Outpatient department, burns, orthopedics and other surgical departments.

Intended population



Maternal / newborn



Specially for baby shower, or swimming, paste on navel



a variety of surgical suture wounds



life / sport careless scratch

WUHAN HUAWEI TECHNOLOGY CO.,LTD.

TRANSPARENT WOUND DRESSING




Specifications



For postoperative wound of orthopedics, slipped disc, extracardial operation and gynaecologic surgery, also used for fixing PICC .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.06.02	10*12	50pcs/box, 20boxes/CTN	40*32*30



For postoperative wound of orthopedics, slipped disc, extracardial operation and gynaecologic surgery, also used for fixing PICC .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.06.06	6*7	100pcs/box, 40boxes/CTN	40*32*30



For postoperative wound of orthopedics, slipped disc, extracardial operation and gynaecologic surgery, also used for fixing PICC .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.06.11	6*7	50pcs/box, 40boxes/CTN	60*22*33



For postoperative wound of orthopedics, slipped disc, extracardial operation and gynaecologic surgery .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.06.12	10*10	25pcs/box, 40boxes/CTN	30*29*27
04.06.13	10*12	25pcs/box, 40boxes/CTN	36*29*27



For postoperative wound of orthopedics, slipped disc, extracardial operation and gynaecologic surgery .


Item No.	Size(cm)	Packing	CTN dimension (cm)
04.06.14	10*15	50pcs/box, 36boxes/CTN	36*38*38
04.06.15	10*20	50pcs/box, 36boxes/CTN	36*38*38
04.06.08	10*25	25pcs/box, 36boxes/CTN	36*38*38
04.06.09	10*30	25pcs/box, 24boxes/CTN	36*38*38



For big ulcer and open wound health care and covering .

Item No.	Size(cm)	Packing	CTN dimension (cm)
04.06.16	15*20	50pcs/box, 20boxes/CTN	36*38*38
04.06.17	15*15	50pcs/box, 20boxes/CTN	36*38*38

- Remarks: 1.All above shape is selectable, we have 4 shapes for each size selection: oval, square, hexagon, U port .
 2.You can select based on your preference. You can also select whether to have absorbent pad or not.
 3.Customized shape and size are welcomed.




WUHAN HUAWEI TECHNOLOGY CO.,LTD.

NONWOVEN WOUND DRESSING



NONWOVEN WOUND DRESSING

Nonwoven wound dressing is made of a soft spunlace non-woven coated with special medical acrylic adhesive as a substrate, the middle of the absorbent pad of cotton fiber has a good moisture permeability, air permeability, good absorb effect, non-adhesive to the wound.



Structure and Characteristics

Ventilation & resistance bacteria: Adopt spunlace non-woven, composite water-based polyurethane film (WPU) non-woven fabrics; cobalt-60 irradiation sterilization, good airtight barrier to bacterial invasion.

Hypoallergenic: advanced coating process, less prone to skin allergies.

Super absorbent: absorbent pads Liquid absorption ≥ 8 times, imported antiadherent film, which guide liquid passing the net type surface, but doesn't stick to the wound. and easy tear off when changing dressing will not cause secondary damage to the wound.

Humanized design: complete specifications, more product styles, we can do according to customers' requirements, provide more choices for clinical use.

Instructions

1. Clean the wound, select the appropriate specification sterile dressing, get it out from the sterile packaging and peel off release paper.
2. The absorbent pad is aligned to wound, gently affixed it to the affected area then thoroughly stripped release paper.
3. Paste the non-woven around smoothly.



Application Postoperative wound, acute and chronic injury wound, small cut wound and scratches, etc.

Clinical departments Outpatient department, burns, orthopedics and other surgical departments.

Intended population



Maternal / newborn

a variety of surgical suture wounds

acute and chronic wounds injury

life / sports careless scratch

Specifications



Applied for incision after care, post-operation wound of laparoscope .

Item No.	Size(cm)	Packing	CTN(cm)
04.10.04	6*7(3*4)	50pcs/box 20boxes/CTN	40.5*21*25.5
04.10.05	7*9		43*25*30.5
04.10.06	9*10(5*6)		43*29*32.5



Applied for postoperative wound of orthopaedic, extracardial and gynaecologic surgery, intervertebral disk, heart surgery .

Item No.	Size(cm)	Packing	CTN(cm)
04.10.11	10*10	50pcs/box 20boxes/CTN	43*29*32.5
04.10.01	10*15	25pcs/box 20boxes/CTN	42.5*29*33
04.10.02	10*20		53*29*33
04.10.13	10*25		57*31.5*32.5
04.10.14	10*30		57*36.5*32.5
04.10.15	10*35		60*37.5*43



Especially applied for urological surgery .

Item No.	Size(cm)	Packing	CTN(cm)
04.10.19	9*20	25pcs/box 20boxes/CTN	53*29*33



Applied for post operative wound of orthopaedic, extracardial and gynaecologic surgery, intervertebral disk, heart surgery .

Item No.	Size(cm)	Packing	CTN(cm)
04.10.07	9*15	25pcs/box 20boxes/CTN	42.5*29*33
04.10.08	9*20		53*29*33
04.10.09	9*25		57*31.5*32.5
04.10.10	9*30		57*36.5*32.5



For eye operatpm after care .

Item No.	Size(cm)	Packing	CTN(cm)
04.10.16	5.7*8.2	50pcs/box 40boxes/CTN	52*29*46
04.10.17	7*9	50pcs/box 20boxes/CTN	



For eye sugery, H shape.

Item No.	Size(cm)	Packing	CTN(cm)
04.10.18	6*12	25pcs/box 20boxes/CTN	39*40*20



For fixing Y or butterfly catheter or catheter with drug kettle .

Item No.	Size(cm)	Packing	CTN(cm)
04.10.21	6*7	50pcs/box 20boxes/CTN	40.5*21*25.5

Remarks: All above shape is selectable our standard size. Customized size and shape are welcomed .

WUHAN HUAWEI TECHNOLOGY CO.,LTD.

Adhesive Bandage and Fixing Series

Catheter Fixing Dressing
Medical Tapes
Silicone Adhesive Tape

CE FDA ISO 9001 ISO 13485



WUHAN HUAWEI TECHNOLOGY CO.,LTD.
CATHETER FIXING DRESSING



CATHETER FIXING DRESSING

The catheter fixing plaster can stabilize the indwelling catheter, reducing the likelihood of traumatic removal, without restricting patient's movement or irritating the skin. all our catheter fixing dressing has special feature unique locking systems that securely anchor the catheter in place, reducing the risk of urethral erosion, bladder spasms and trauma .

Specifications



Nasal feeding tube fixation device

Item No.	Size(cm)	Packing
09.03.01	2*7	50pcs/box



PICC catheter fixing device

Item No.	Size(cm)	Packing
09.03.02	2*8	50pcs/box



I.V. Fixing device

Item No.	Size(cm)	Packing
09.03.03	3*9	50pcs/box



PICC catheter fixing device

Item No.	Size(cm)	Packing
09.03.04	3*9	50pcs/box




PICC catheter fixing device

Item No.	Size(cm)	Packing
09.03.05	3*9	50pcs/box



Epidural fixation device

Item No.	Size(cm)	Packing
09.03.07	3*9	50pcs/box




I.V. Fixing device

Item No.	Size(cm)	Packing
09.03.08	6*7	50pcs/box



urethral catheter fixing device

Item No.	Size(cm)	Packing
09.03.09	3*9	50pcs/box



urethral catheter fixing device

Item No.	Size(cm)	Packing
09.03.10	8*12	50pcs/box



urethral catheter fixing device

Item No.	Size(cm)	Packing
09.03.11	6*11	50pcs/box



urethral catheter fixing device

Item No.	Size(cm)	Packing
09.03.12	6*11	50pcs/box




nasogastric feeding tube, Fixing device

Item No.	Size(cm)	Packing
09.03.14	3*8	100pcs/box



nasogastric feeding tube, Fixing device

Item No.	Size(cm)	Packing
09.03.15	7*10	100pcs/box



Tracheal intubation fixation

Item No.	Size(cm)	Packing
09.03.16	12*8	50pcs/box

MEDICAL TAPE

Medical tape is heat-seal coated by the pressure sensitive adhesive on the non-woven, release paper, PE film, cotton cloth, polyester cotton cloth, polyester cloth, artificial silk, silk, imitation silk, imitation cotton and other raw materials . Solid paste, not easy to fall, no harm to skin, no adhesive remains to skin when peeling off .

Cotton tape Non-woven tape PE tape Silk Tape



Tape Name	Features/Advantages
Cotton tape (Plaster)	Strong viscosity, thick and difficult to tear.
Non-woven tape (Micro-porous and breathable)	Hypoallergenic, breathable and easy to tear.
PE tape (Transparent)	Transparent and breathable, easy to tear and good viscosity
Silk Tape	Super hypoallergenic silk, good viscosity

Product Features

Strong viscosity, easy to tear: good viscosity, there is no pain to skin when peeling off .

Hypoallergenic: hypoallergenic formula, no harm to skin .

Porous breathable: good breathability, soft and thin just like a second skin of human .

Wide application: for fixing plaster, gauze, needles, sleeves, etc .

PE TAPE



Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.05.15	1.25*914	24rolls/box 30boxes/CTN	35*34.5*28.5	Used in surgical drapes fixation, infusion plaster fixation and banding and fixation of other medical device.
04.05.16	2.5*914	12rolls/box 30boxes/CTN		
04.05.17	5.0*914	6rolls/box 30boxes/CTN		
04.05.18	7.5*914	4rolls/box 30boxes/CTN		

PAPER TAPE



Item No.	Size(cm)	Packing	CTN (cm)	Sterile Method
04.05.19	1.25*914	24rolls/box 30boxes/CTN	28.5*28*24	Used in surgical drapes fixation, infusion plaster fixation and banding and fixation of other medical device.
04.05.20	2.5*914	12rolls/box 30boxes/CTN		
04.05.21	5.0*914	6rolls/box 30boxes/CTN		
04.05.22	7.5*914	4rolls/box 30boxes/CTN		

WUHAN HUAWEI TECHNOLOGY CO.,LTD.

MEDICAL TAPE



SILK TAPE



Item No.	Size	Packing	CTN (cm)	N/W (kg)
04.05.23	1.25cm*5m	12pcs/box 60boxes/CTN	56*25*29	14
04.05.24	1.5cm*5m			14
04.05.25	2.5cm*5m	12pcs/box 30boxes/CTN		13
04.05.26	5cm*5m	6pcs/box 30boxes/CTN		13
04.05.27	7.5cm*5m	6pcs/box 24boxes/CTN	51.5*24*32.5	19
04.05.28	10cm*5m	6pcs/box 24boxes/CTN	46*34.5*32.5	20
04.05.29	1.25cm*10yds	24pcs/box 48boxes/CTN	45*32.5*33.5	19
04.05.30	2.5cm*10yds	12pcs/box 48boxes/CTN		19
04.05.31	5cm*10yds	6pcs/box 48boxes/CTN		19
04.05.32	7.5cm*10yds	4pcs/box 48boxes/CTN		19

ZINC OXIDE ADHESIVE PLASTER(SIMPLE BOX PACKAGE)



Item No.	Size	Packing	CTN (cm)	N/W (kg)
04.05.35	1.25cm*5m	24pcs/box 30boxes/CTN	56*25*29	14
04.05.36	2.5cm*5m	12pcs/box 30boxes/CTN		13
04.05.37	5cm*5m	6pcs/box 30boxes/CTN		13
04.05.38	7.5cm*5m	6pcs/box 24boxes/CTN	51.5*24*32.5	19
04.05.39	10cm*5m	6pcs/box 24boxes/CTN	46*34.5*32.5	20

ZINC OXIDE ADHESIVE PLASTER(PLASTIC TUBE PACKAGE)



Item No.	Size	Packing	CTN (cm)	N/W (kg)
04.05.40	1.25cm*5m	24pcs/box 30boxes/CTN	37*36.5*42.5	17
04.05.41	2.5cm*5m	12pcs/box 30boxes/CTN		15
04.05.42	5cm*5m	6pcs/box 30boxes/CTN	37*36.5*42.6	15
04.05.43	7.5cm*5m	6pcs/box 24boxes/CTN		19
04.05.44	10cm*5m	4pcs/box 30boxes/CTN	46*34.5*32.5	24

Remarks: 1. Above are all standard size, for customized size and shape, material, color and package are available upon request.
2. Each size can be made into different material, the picture is just for your reference of size and shape.

SILICONE ADHESIVE TAPE

It is made of PU film or nonwoven coated with silicone gel. It is applicable to variety of in-dwelling catheters, electrodes, or other external fixation that requires gentle but robust skin protection and adhesion. Protect the skin which have the risk of ulcer of catheters, or the skin need to be radiation therapy. The tape can be torn into the size you need, as there is cutting line on the film (or tearable nonwoven).


Advantages

1. Silicone tape make an excellent addition to a medical emergency supply kit .
2. This Tape is designed to not harm the skin when it is being removed .
3. This Tape is adhesive enough that it can be worn in the shower without losing effectiveness .
4. single use silicone tape , Waterproof silicone tape, but permeable .


Specifications



Item No.	Product Name	Size	Package	CTN(cm)
11.05.01	Silicone tape deluxe type	2.5cm*1.5m	1roll/box, 200boxes/CTN	38*33*29
11.05.02		3.5cm*1.5m	1roll/box, 200boxes/CTN	38*33*29
11.05.03		4cm*1.5m	1roll/box, 200boxes/CTN	38*33*29
11.05.04		2.5cm*5m	1roll/box, 200boxes/CTN	38*43*39
11.05.05		3.5cm*5m	1roll/box, 200boxes/CTN	38*43*39
11.05.06		4cm*5m	1roll/box, 200boxes/CTN	38*43*39




Item No.	Product Name	Size	Package	CTN(cm)
11.06.01	Silicone tape waterproof- type	2.5cm*1.5m	12rolls/box, 48boxes/CTN	47*42.5*32
11.06.02		3.5cm*1.5m	12rolls/box, 48boxes/CTN	47*42.5*32
11.06.03		4cm*1.5m	12rolls/box, 48boxes/CTN	47*42.5*32
11.06.04		2.5cm*5m	12rolls/box, 48boxes/CTN	47*42.5*32
11.06.05		3.5cm*5m	12rolls/box, 48boxes/CTN	47*42.5*32
11.06.06		4cm*5m	12rolls/box, 48boxes/CTN	47*42.5*32



Item No.	Product Name	Size	Package	CTN(cm)
11.07.01	Silicone tape Economic- type	2.5cm*1.5m	1roll/box, 400boxes/CTN	33.5*32*22
11.07.02		3.5cm*1.5m	1roll/box, 400boxes/CTN	42*33.5*22
11.07.03		4cm*1.5m	1roll/box, 400boxes/CTN	52*33.5*22
11.07.04		2.5cm*5m	1roll/box, 400boxes/CTN	49*31*32
11.07.05		3.5cm*5m	1roll/box, 400boxes/CTN	49*40*32
11.07.06		4cm*5m	1roll/box, 400boxes/CTN	49*48*32

WOUND SKIN CLOSURE STRIP

Wound closure strip are made of non-woven, rayon-backed material reinforced with filaments for strength coated with Hypoallergenic adhesive. Provides wound support and assists in increasing the tensile strength of the wound .



Application

1. For closure of general surgery wounds .
2. For closure of gynecologic surgery wound such as caesarean section surgery wound .
3. For wound care after immediate suture or stitch removal .

Usage

1. Clean the wound with saline water and dry the wound area .
2. Stick the wound skin closure to the wound vertically from the middle of the wound .
3. Apply the wound skin closure from the two ends of the wound to the middle part to avoid wound opening .
4. Stick to the rest parts of wound .
5. The medical band aid can be left in place for several days and is sticking well .




Specifications

Item No.	Size(mm)	Packing	Inner Pouch Size(cm)	Inner Box size(cm)	CTN(cm)
03.01.01	3*75	5pcs/sheet, 1sheet/bag, 100bags/box, 40boxes/CTN	6.3*14.2	14.7*5.8*6.8	31.5*31*29.5
03.01.02	6*38	6pcs/sheet, 1sheet/bag, 100bags/box, 40boxes/CTN	6.3*14.2	14.7*5.8*6.8	31.5*31*29.5
03.01.03	6*75	3pcs/sheet, 1sheet/bag, 100bags/box, 40boxes/CTN	6.3*14.2	14.7*5.8*6.8	31.5*31*29.5
03.01.04	6*100	5pcs/sheet, 2sheets/bag, 50bags/box, 30boxes/CTN	7.3*16.6	17*5.8*7.7	36*31*25
03.01.05	12*50	3pcs/sheet, 2sheets/bag, 50bags/box, 40boxes/CTN	6.3*14.2	14.7*5.8*6.8	31.5*31*29.5
03.01.06	12*100	3pcs/sheet, 2sheets/bag, 50bags/box, 30boxes/CTN	7.3*16.6	17*5.8*7.7	36*31*25
03.01.07	12*125	3pcs/sheet, 2sheets/bag, 50bags/box, 30boxes/CTN	7.3*17.5	18*5.8*7.7	38*31*25
03.01.08	25*100	2pcs/sheet, 2sheets/bag, 50bags/box, 30boxes/CTN	9*16.6	17*5.8*9.5	36*31*30.5

WUHAN HUAWEI TECHNOLOGY CO.,LTD.

WOUND SKIN CLOSURE STRIP



WUHAN HUAWEI TECHNOLOGY CO.,LTD.

Personal Healthcare Series

Silicone Scar Gel Sheet

Acne Patch

Pain Relief Patches


Cold Therapy Products

Functional Herbal Patches

CE FDA ISO 9001 ISO 13485

WUHAN HUAWEI TECHNOLOGY CO.,LTD.

SILICONE SCAR GEL SHEET



SILICONE SCAR GEL SHEET

Silicone Scar Gel Sheet has similar effects to the human skin, with the character of elasticity, moisture permeability, air permeability, pressure sensitive resistance, water repellency and has barrier function of bacteria. The heterogeneous collagen layer which contacts to scar tissue is rich in collagen and with good biological compatibility, can stick on the scar tissue, and no other dressings are needed to be fixed .

Main Principle

The scar extrudes from the normal skin in irregular shape and by the form of red and hard benign lumps and even with contracture, itching ordains, In case of excessive growth, the scar might affect the appearance or lead to dysfunction or even cancer .


Performance

Silicone gel can directly act on the scar surface. As a semi-permeable membrane, it can partly restrict the water permeability on the scar surface and promote the moisture accumulation under the scar epidemics, resulting in hydration. Studies have shown that using the Silicone gel patch can greatly reduce the number of the positive cells of the fibroblast1 growth factor (TGF 1), which indicates that the Silicone gel patch acts on scars possibly by adjusting the relationship between cell factors and fibroblasts. Experiments show that the silicone oil slowly released from the Silicone gel membrane can inhibit the strengthening of the fibroblasts inside the scar tissue, thereby inhibiting the formation of hypertrophic scars. Through the trilogy of "softening, restoring and regeneration", it can awaken the memory of the skin to fully repair various scar tissues .

Functions

1. Lock moisture and pressure on the skin
2. Scar Remove
3. Reduce scar formation

Specifications



Item No.	Size(cm)	Package	CTN(cm)	Remarks
11.08.01	2*12	1piece/box, 200boxes/CTN	52*39*36	The color of the sheet is changable.
	4*12	1piece/box, 200boxes/CTN	52*39*36	
	6*12	1piece/box, 200boxes/CTN	52*39*36	



WUHAN HUAWEI TECHNOLOGY CO.,LTD.

ACNE PATCH

ACNE PATCH

Acne patch visibly shows it's working by absorbing pus and fluids , Drug-free and without harsh chemicals or drying formulas, it is cruelty-free and vegan friendly ,with advanced hydrocolloid material, fast absorption on exudates from acnes, it protects blemishes for faster healing overnight, pimple patches act as a protective cover and prevent picking or touching, effectively reduce the redness of acnes, and help in offering a moist environment, accelerate the acne healing without scars left .

Application

1. Especially adaptable to broken acne with exudates .
2. For shallow and small wounds caused by injury .
3. For wounds caused by laser and other minor surgery .
4. For protection of small wounds on face, avoiding the blockage and infection caused by make-up .


Specifications



Item No.	Specification	Sheet(cm)	Packing	Remarks
04.19.01	Dia14mm*16 dots	7*7	1sheet/pouch, 10 pouches/box, or mixed size packing	0.06mm thickness, ultra thin
04.19.02	Dia12mm*16 dots	6.5*6.5		
04.19.03	Dia10mm*16 dots	6*6		
04.19.04	Dia8mm*16 dots	5*5		
04.19.05	4.5*4.5cm*2 dots	5*5		




Item No.	Specification	Sheet(cm)	Packing	Remarks
04.19.06	Dia12mm*12 dots	6*7.5	1sheet/pouch, 2-10 pouch/box, or mixed size packing	border 0.4mm, middle 0.6mm
04.19.07	Dia12mm*12 dots	6*7	1sheet/pouch, 2-10 pouch/box, or mixed size packing	0.5mm thickness,
04.19.08	Dia12mm*12 dots	6*7		
04.19.09	Dia12mm*12 dots	6*7		



Item No.	Specification	Sheet(cm)	Packing	Remarks
04.19.10	Dia12mm*8 dots Dia 10mm*8dots Dia 8mm*8dots	7*8	1sheet/pouch, 2-10 pouch/box, or mixed size packing	0.5mm thickness,
04.19.11	Dia10mm*12 dots Dia 8mm*24 dots	6.5*9		



WUHAN HUAWEI TECHNOLOGY CO.,LTD.
PAIN RELIEF PATCHES



PAIN RELIEF PATCHES


Structure and Performance

Pain is a worldwide problem. Everyone will face all kinds of pain from birth to death. Our company provides all kinds of pain solutions. We are the largest production base of plaster in China. Here are some of our most conventional pain-killing patches :


Specifications



Item	Material /structure	Applicable scope	Speicification of item(mm)	Packaging
Fuyang moxibustion patch	Far Infrared magnetic powder, medical adhesive, non-woven fabric, menthol, ginger juice,	Suitable for lumbar disc herniation, ankylosing spondylitis, asthma, frequent urination, irregular menstruation, Pain and discomfort caused by driving and poor health.	550*80	1 pc/pouch, 3 pouches/box, 80 boxes/CTN



Item	Material /structure	Applicable scope	Speicification of item(mm)	Packaging
Infrared pain relief patch	Far Infrared magnetic powder, medical adhesive, non-woven fabric, menthol	Suitable for pains of muscle and joint resulting from the Wind-Cold-Dampness (Chinese-called pathogenic factors), and the pains of soft-tissue contusion.	100*70	2 pcs/pouch, 3 pouches/box, 120 boxes/CTN



Item	Material /structure	Applicable scope	Speicification of item(mm)	Packaging
Traditional black magnet plaster for pain relief	Euphorbia herb, gallnut, wild chrysanthemum flower, camphor, permanent magnet	Applies to humans accupoints as anti-inflammatory and pain relief, use for bruises, muscle strain, arthritis, adjuvant joint pain and numbness, frozen shoulder, osteoarthritis, aching pain and other pain disorders subsidiary treatment	120*160	1pc/pouch, 4 pouches/box, 90 boxes/CTN



Item	Material /structure	Applicable scope	Speicification of item(mm)	Packaging
Capsicum plaster	Capsicum: Can expanding the capillaries of skin, thus to promote blood circulation, has a remarkable effect on closed soft tissue injury. Menthol: Menthol can restrain tumefaction of rat in the test., Relieve local pain rapidly and powerful. Camphor: also has a good effect on analgesic and anti-inflammatory.	Uitable for relief of all types of minor aches; include muscle pain, shoulder and back pain, joint pain, strain, bruises and sprains, etc	100*70	1 pc/pouch, 24 pouches/box, 180 boxes/CTN
			120*80	1 pc/pouch, 24 pouches/box, 180 boxes/CTN
			160*100	1 pc/pouch, 24 pouches/box, 120 boxes/CTN
			160*120	1 pc/pouch, 24 pouches/box, 120 boxes/CTN
			180*120	1 pc/pouch, 24 pouches/box, 80 boxes/CTN



Item	Material /structure	Applicable scope	Speicification of item(mm)	Packaging
Analgesic fomentation plaster	Far Infrared magnetic powder, medical adhesive, non-woven fabric, heating bag, herbal patches	Suitable for Rheumatic and algidity, cervical spondylosis, scapulohumeral periarthritis, lumbar muscle strain, osteoarthritis of bone hyperplasia, sciatica and injuries of the soft tissue.	188*88	1 pc/pouch, 3 pouches/box, 80 boxes/CTN




Item	Material /structure	Applicable scope	Speicification of item(mm)	Packaging
Self heating pad warmer	Iron, vermicolite, salt, activated carbon, medical adhesive, non-woven fabric, paper	1. Effectively soothe pains from rheumatism, muscular aches and pains, joint stiffness, backache, strains and sprains, pain caused by cold Stomachache. 2. Provide menstrual care for female and abate pains form dysmenorrheal (menstrual cramps), endometriosis. 3. Can be used by athletes, outdoor enthusiasts, sportsmen, spectators, skiers and construction worker anyone who hate cold weather.	133*100	1pc/pouch, 3 pouches/box, 100 boxes/CTN

COLD THERAPY PRODUCTS


Structure and Performance

Hydrogel is composed of water as dispersion medium and some materials such as polymer and crosslinking agent. It is good skin-friendly, gentle and moderate, does not hurt the skin, widely range of application, according to the different parts of the application, with different essential oil composition in the patch, it can achieve different effects. Hydrogels are rich in water, and the evaporation of water will take away a lot of heat, thus playing a cooling effect and effectively reducing the local temperature of the skin. Cold compress can make capillary systole, reduce local congestion, can make the sensitivity of nerve endings reduces and reduce pain, cool down and act as an antipyretic, can reduce local blood flow, prevent inflammation and suppurative diffusion. It can dissipate heat conduction in the body, increase heat dissipation and lower body temperature. At present, hydro gels are widely used in the medical field. In terms of personal care, we have follow products for our customer choice, we also can do our customer's formulation and add any ingredients if they like to achieve special treating purpose .


Specifications



Item	Material /structure	Applicable scope	Speicification of item(cm)	Packaging
Hydrogel pain relief patch	Backing layer, gel layer (contains medical gel, glycerol, water, menthol, camphol, capsicum) and excellent release layer	Cooling relief from bruises, inflammation, sprains and strains, Chronic joint or muscle pain, easy to Keep it in your car, sports bag, first aid kit, or anywhere you may need access to instant cooling support.	12*8 10*14 10*20	2 pcs/pouch, 2 pouches/box, 120 boxes/CTN



Item	Material /structure	Applicable scope	Speicification of item(cm)	Packaging
Hot and cold patch	Backing layer, gel layer (contains medical gel, glycerol, water, menthol) and excellent release layer	Temporarily relieves minor aches and pains of muscles and joints due to arthritis, simple backache, sprains, bursitis, tendonitis, bruises, cramps	7*10 8*12 10*14 10*20	2 pcs/pouch, 2 pouches/box, 120 boxes/CTN



Item	Material /structure	Applicable scope	Speicification of item(cm)	Packaging
Cooling gel patch	Backing paper, medical gel, glycerol, water, menthol, lavender oil, Eucalyptus Oil, and excellent release paper	It is applicable to physically reduce the local temperature and prevent fever heat, heat convulsions to protect the brain effectively. It can be used to relieve local pain caused by trauma without damage. And it also can be used to reduce the local temperature of acute febrile, headache, toothache, traumatic injury and isolation.	10*4 11*4.5 12*5 10*14	1 pc/pouch, 4 pouches/box, 240 boxes/CTN

WUHAN HUAWEI TECHNOLOGY CO.,LTD.

COLD THERAPY PRODUCTS



FUNCTIONAL HERBAL PATCHES

We are the largest manufacturer of gluing and adhesive products in China. We can add any essential oil in the glue to achieve special functions and curative effects. Here are some personal care products we have produced at present .


Specifications




Item	Material /structure	Applicable scope	Size(mm)	Packing
Infrared cough relief patch for adult	Far-infrared materials, medical adhesive, non-woven fabric herbal extracts	Melt and eliminate phlegm, relieve cough, comfort throat, clear lung, and improve immunity, special for adult	80*55 (38*30)	4 pcs/pouch, 2 pouches/box, 120 boxes/CTN
			60*50 (38*30)	



Item	Material /structure	Applicable scope	Size(mm)	Packing
Pharyngitis patch	Far-infrared materials, medical adhesive, non-woven fabric herbal extracts	Melt the phlegm and relieve cough, relieve pharyngitis, comfort throat, clear lung, and improve immunity for health protection.	40*50	2 pcs/pouch, 2 pouches/box, 240 boxes/CTN



Item	Material /structure	Applicable scope	Size(mm)	Packing
Nose Pore patch	Deionized water, carbon powder, poval, hamamelis virginiana extract, hydroxymethyl cellulose, glycolic acid, allantoin, essence	Remove blackhead, makeup residues of harmful substances, and strengthen the effect to absorb blackhead. deep clean and nourish skin both inside and outside, promote the skin metabolism, active skin cells, all of the functions lead to make the skin white and luster.	39*90	1 pc/pouch, 10 pouches/box, 40 boxes/CTN



Item	Material /structure	Applicable scope	Size(mm)	Packing
Rhinitis patch	Semen brassicae, euphorbia kansui, manchurian wildginger, lily magnolia, syzygium aromaticum, rhizoma corydalis, cocklebur, rhizoma zingiberis, cinnamon and mucilage	Used to for simple rhinitis, hypertrophic rhinitis, allergic rhinitis, chronic rhinitis, sinusitis rehabilitation, also used to cure and recover nasal congestion caused by rhinitis and colds, sneezing, nasal itching, excessive nasal mucus, Nose turbidity, nasal swelling, hyposmia, nasal heavy, headache, dizziness, red eye, tearing, and sleep disorders, etc	102*30 90*44	2 pcs/pouch, 5 pouches/box, 240 boxes/CTN


WUHAN HUAWEI TECHNOLOGY CO.,LTD.

FUNCTIONAL HERBAL PATCHES





Item	Material /structure	Applicable scope	Size(mm)	Packing
Foot Detoxification Patch	Wooden vineger powder, bamboo vineger powder, tourmaline, chitin, vitamin C, plant powder, onion powder, cornstarch and other ingredients	1.Dispelling toxins accumulated in absorbent system, promoting functions of absorbent system and strengthening the immunity, 2. Promoting blood circulation & metabolism,3. Activating cells, improving functions of vital organs, 4. Relieving foot fatigue,5. Relaxing muscles & tendons and eliminating moisture, replenishing vital essence and strengthening kidney,6. Relieving swelling & pain, 7. Promoting sleeping	4.5g ±0.5g/ patch, 20 patches/ box	1pc/pouch, 10 pouches/box, 120 boxes/CTN



Item	Material /structure	Applicable scope	Size(mm)	Packing
Body slim patch	Pinellia, lotus leaf, tuckahoe, alisma, hawthorn, focal three fairy, cassia seed, patchouli essential oil	Use the medicine to stimulate the acupoint, regulating vital energy and activating blood, eliminating stasis in the body, recovering normal metabolism function in body, adjust gastrointestinal tract and pancreas function, conditioning internal secretion, activate lipase activity speed up lipolysis, and further resolve and consume excessive storing stasis and fat, and make it expel from the body accompany with the metabolism.	70*70	1 pc/pouch, 30 pouches/box, 20 boxes/CTN



Item	Material /structure	Applicable scope	Size(mm)	Packing
Corn plaster	Salicylic acid, safflower, borneol, Sessile Stemona Root, Java Brucea Fruit	Relieves pain immediately & helps remove the corn	65*22	6 pcs/box, 500 boxes/CTN



Item	Material /structure	Applicable scope	Size(mm)	Packing
Mosquito patch	Laminated non-woven, Citronella oil, lemon Oil, Eucalyptus oil, protection paper	Respell the mosquito around the wearer, also with decoration effect	Whole sheet 74*110	6 pcs/sheet, 1 pc/pouch, 6 pouches/box
			Patch size 34*34	



Item	Material /structure	Applicable scope	Size(mm)	Packing
Bian stone body slim patch	Bian stone	1. Conditioning organs, spleen and stomach function, so that it tends to normal, regulate the pathological appetite, reduce excess calorie intake, 2. Increase the body's metabolic function, dispelling phlegm and eliminating dampness, 3. Bian stone needle slimming patch through the physical effects of the composition of the stone to stimulate the meridian points to stimulate the innate meridian system to enhance the function of the body organs	100*100	1 sheet/ pouch, 14 pouches/ box



Item	Material /structure	Applicable scope	Size(mm)	Packing
Hangover patch	Peppermint, borneol, Piper methysticum extract, green tea extract, taurine and so on	1. Improve concentration, ease the brain fatigue, prevent excessive using brain and improve the learning efficiency 2. who use computer long time: get his brain refresh 3. who sit long time, such as driver: relieve fatigue, can also be used for the drunken and unwell, and so on.	30*30	6 pcs/ sheet, 1 sheet/ pouch, 5 pouches/ box



Item	Material /structure	Applicable scope	Size(mm)	Packing
Stop smoking patch	Codonopsis pilosula, chrysanthemum, pepper mint, honeysuckle, clove, Natural tobacco	1. Relief of withdrawal symptoms associated with giving up smoking. 2. Reduce the frequency of smoking by using the patches daily. 3. Reduce the cravings for smoking.	50*50	2 pcs/ pouch, 6 pouches/ box

Remarks

1. OEM and ODM are welcome .
2. The above patches can be added with extra ingredients according to your requirements, or use our standard ingredients .



Isolation Gowns

Isolation gowns are examples of personal protective equipment used in health care settings. They are used to protect the wearer from the spread of infection or illness if the wearer comes in contact with potentially infectious liquid and solid material. They may also be used to help prevent the gown wearer from transferring microorganisms that could harm vulnerable patients, such as those with weakened immune systems. Isolation gowns are one part of an overall infection-control strategy.



Item No.	Item Name	Specification	Packing	CTN (CM)	G.W.	Feature
20.03.01	Isolation gown	(S,M,L,XL,XXL)	10pcs/bag, 100pcs/ctn	52*32*50	7.5KG	PP+PE, PP,SMS

Uses of isolation gown: Operating room, clinical application, first aid room



- ✓ Features of isolation gown:
- ✓ Flexible soft material, which is very comfortable to wear
- ✓ Economic, light weight, comfortable, breathable
- ✓ Sleeves with knitted cuff or elastic cuff
- ✓ Moisture repellent
- ✓ Latex free

Application of isolation gown: Operating room, clinical application, first aid room etc.

WUHAN HUAWEI TECHNOLOGY CO.,LTD.

PERSONAL PROTECTIVE EQUIPMENT




FACE SHIELD



- ✓ Medical-grade protective material
- ✓ 180° facial isolation protection
- ✓ Blocking droplets, aerosols
- ✓ The panel is light and detachable

Product features

1. The panel is made of PET material, colorless and transparent, anti-fog, and can effectively protect the face skin and eyes
2. It has the characteristics of double-sided anti-fog, anti-static and anti-vertigo, which can effectively prevent the blurred vision caused by temperature difference and water vapor and has no stimulating effect on the eyes.
3. The panel is thin and light in weight, comfortable to wear.
4. Effectively block the potential patient's blood, body fluids and harmful chemical spatter and other harm to human body



Item No.	Item Name	Specification	Packing	CTN (CM)	G.W.	Feature
20.02.01	face shield	33*22cm	1pcs/bag, 200pcs/ctn	35*49*45cm	9.3kg	transparent

Application scope


Use widely and meet the daily need

It is used in health protection; resistance of coughing and sneezing, dust particles; Riding in rainy days

Disposable Face Mask

Main components

The mask consists of 3 Ply layered filter (Lining: Nonwovens; Filter layer: Melt-blown fabric, outer layer: Nonwovens), nose clip and mask belt.



BFE≥95% , Soft , Breathable

Scope of application

Suitable for covering the user's mouth, nose and jaw, and providing the wearer's health care in a general environment. Create a physical barrier to guard particles, smog, pollen.

Method of application

1. Unfold the folds to make the mask fully extended; put on your face and pressed the nose clip close to the bridge of your nose.
2. Cover your face while wearing an ear loop; Pull the mask down to the chin, making sure a large area of protection;
3. Always check when you wear the mask: put your hands over your mask, then exhale. If you feel air leaking from the nose clip, you need to re-adjust it.
4. The best effective time of the mask is 4 hours.



Item No.	Item Name	Specification	Packing	CTN (CM)	G.W.	Feature
20.01.01	Disposable Face Mask	17.5*9.5cm	50pcs/box, 1100pcs/ctn	56*30*35cm	5kg	BFE≥95%

Matters need attention

1. Disposable product, cannot produce oxygen and cannot be used in anoxic environment.
2. This product is strictly prohibited to be used in the environment contain gas or harmful dust.
3. This product is normal grade and Non-sterile. If the package is damaged, please do not use it.
4. Please keep away from fire during use.

WUHAN HUAWEI TECHNOLOGY CO.,LTD.

Medical Raw Material Jumbo Roll Series

PU Material


PE Material

Nonwoven and Woven Material

Foam Material

Hydrocolloid Dressing Roll

CE FDA ISO 9001 ISO 13485



WUHAN HUAWEI TECHNOLOGY CO.,LTD.

MEDICAL COATING RAW MATERIAL SERIES



MEDICAL COATING RAW MATERIAL SERIES

We supply many kinds of raw or semi-finished material such as PU, PE, elastic fabric, cotton fabric or non-woven coating material, meaning while, we can supply the material based on customers' requirement.

With our help, you can set up the fast standard production line for all medical products.

All of products are in medical grade, our medical coated raw materials has following character:

1. Various material, For example:

Non-woven: we have follow types with different weight and colors for coating: Spun-laced non-woven, thermal-bond nonwoven, melt-blown pp nonwoven, Hot-blast non woven.









PU film: With PET backing, With PE backing, Without any backing.



PE film: Frost opaque outlooks, transparent smooth out looks, with micro-perforation permeability, With color imprinted designs.

2. Various coating technology: strip-coated, grid-coated.











3. Unique adhesive and wide choice of adhesives type: such as acrylic adhesive, hot medical adhesive, hydrogel adhesive, silicone adhesive, etc.

PU Material

	Item name	Fabric weight/ thickness	Size/roll	Packing	
	PU film(for wound dressing, I.V dressing, etc) with S cutting liner	0.04mm PET +0.02mm PU + adhesive + 80g/m2 protection paper	10cm*300m (any width and length ok)	20 rolls/CTN	
	Item name	Fabric weight/ thickness	Size/roll	Packing	
	PU film(for wound dressing, I.V dressing, etc) with cutting liner	0.05mm PE + 0.02mm PU + adhesive + 80g/m2 protection paper	10cm*300m (any width and length ok)	20rolls/CTN	
	Item name	Fabric weight/ thickness	Size/roll	Packing	
	PU film with finger lift(for surgical incisive dressing, wound dressing, etc)	0.05mm PE + 0.02mm PU + adhesive + 80g/m2 protection paper	10cm*300m (any width and length ok)	20 rolls/CTN	
	Item name	Fabric weight/ thickness	Size/roll	Packing	
	PU film with top paper backing	110g CCK + 0.02mm PU + adhesive + 80g/m2 protection paper	10cm*300m (any width and length ok)	20 rolls/CTN	

	Item name	Fabric weight/ thickness	Size/roll	Packing	
	Iodophor PU film(for surgical incise dressing)	0.02mm PU + adhesive with iodophor + 80g/m2 protection paper	45cm*300m (any width <48cm)	1 roll/CTN	


PE Material

	Item name	Fabric weight/ thickness	Size/roll	Packing	
	PE transparent film(can be used for band aid, also for surgical dressing)	0.02mm PE + adhesive + 80g/m2 protection paper	50cm*300m (any width and length ok)	20 rolls/CTN	
	Item name	Fabric weight/ thickness	Size/roll	Packing	
	PE transparent film with finger lift (for surgical incisive dressing)	0.02mm PE + adhesive + 80g/m2 protection paper	20cm*300m (any width and length ok)	20 rolls/CTN	
	Item name	Fabric weight/ thickness	Size/roll	Packing	
	PE tape	72g/square meter	54cm*500m (any width and length ok)	1 roll/CTN	
	Item name	Fabric weight/ thickness	Size/roll	Packing	
	Transparent high elastic PE tap (For band aid)	0.06mm thick	7.6cm*500m (any width and length ok)	20 rolls/CTN	
	Item name	Fabric weight/ thickness	Size/roll	Packing	
	Skin color micro-perforated PE material (for band aid)	0.08mm thick	7.6cm*500m	1 roll/CTN	




Item name	Fabric weight/thickness	Size/roll	Packing
Net type PE tape (for band aid, bottle mouth protection patch)	0.06mm thick	7.6cm*500m	20 rolls/CTN







Item name	Fabric weight/thickness	Size/roll	Packing
Iodophor PE film (for surgical incise dressing)	0.02mm PE + adhesive with iodophor + 80g/m ² protection paper	45cm*300m (any width <48cm)	1 roll/CTN




Nonwoven and Woven Material





Item name	Fabric weight/thickness	Size/roll	Packing
Non-woven roll (For fixing up, bandage or wound dressing)	45g/square meter	50cm*300m (any width and length ok)	1 roll/CTN







Item name	Fabric weight/thickness	Size/roll	Packing
Zinc oxide white cotton (For zinc oxide medical tape)	115g/m ²	91cm*50m	1 roll/CTN







Item name	Fabric weight/thickness	Size/roll	Packing
Silk tape	80g/m ²	50cm*300m	1 roll/CTN







Item name	Fabric weight/thickness	Size/roll	Packing
Paper tape	65g/m ²	54cm*500m	1 roll/CTN






Item name	Fabric weight/thickness	Size/roll	Packing
Skin color elastic fabric material (for band aid)	115g/m ²	7.4cm*300m (any width and length ok)	20 rolls/CTN







Item name	Fabric weight/thickness	Size/roll	Packing
Capsicum plaster jumbo roll	Perforated cotton + zinc oxide adhesive with herbals oils + protection paper	12cm*300m	4 rolls/CTN




Hydrocolloid Dressing Roll




Item name	Fabric weight/thickness	Size/roll	Packing
Hydrocolloid dressing roll	0.4mm, 0.5mm, 0.6mm, 0.8mm, 1.0mm, 1.2mm, Thickness hydrocolloid + protection paper	23cm*100m	20 rolls/CTN



Foam Material



Item name	Fabric weight/thickness	Size/roll	Packing
Foam dressing no adhesive	3mm thickness Foam laminated with 0.02mm PU	28cm*100m	1 roll/CTN
Foam dressing no adhesive	5mm thickness laminated with 0.02mm PU	28cm*100m	1 roll/CTN



Remarks: 1. We can do adhesive coating on all type materials .
2. All of products can be any width and length depend on your requirements .